

---

# Signal Server – UserManual

*Release 1.0*

**Christian Breitwieser – c.breitwieser@tugraz.at**

January 12, 2011

## Contents

<b>1</b>	<b>Introduction</b>	<b>ii</b>
1.1	Supported hardware . . . . .	ii
1.2	Planned hardware support . . . . .	ii
1.3	Planned extensions . . . . .	ii
1.4	License . . . . .	ii
1.5	Contact . . . . .	ii
<b>2</b>	<b>Hardware Requirements</b>	<b>iii</b>
<b>3</b>	<b>Installation Instruction</b>	<b>iii</b>
3.1	Hardware Drivers . . . . .	iii
3.2	Debian/Ubuntu . . . . .	iii
3.3	Windows . . . . .	iii
<b>4</b>	<b>Using the Signal Server</b>	<b>iii</b>
4.1	Generic Usage . . . . .	iii
4.2	Debian/Ubuntu . . . . .	iv
4.3	Windows . . . . .	iv
<b>5</b>	<b>The XML Configuration File – Basics</b>	<b>iv</b>
5.1	Subject Information . . . . .	iv
5.2	Server Settings . . . . .	iv
5.3	The Hardware sections . . . . .	iv
	Mode . . . . .	iv
	Device Settings . . . . .	iv
	Channel Settings . . . . .	v
<b>6</b>	<b>Specific Hardware Configuration with the XML Config File</b>	<b>v</b>
6.1	Sine Generator . . . . .	v
6.2	g.USBamp . . . . .	v
6.3	g.Mobilab . . . . .	v
6.4	g.BSBamp . . . . .	v
6.5	Brainamp Series . . . . .	v
6.6	EEG Simulator . . . . .	v
6.7	Generic Joysticks . . . . .	v
6.8	IntegraMouse + Generic Mouses . . . . .	v

---

# 1 Introduction

The TOBI signal server describes a program using TiA (TOBI interface A) to acquire and distribute raw biosignals. Implementation was done using C++, static and shared library packages are available for Debian based systems (32bit and 64 bit), Microsoft Windows Xp and Windows 7.

## 1.1 Supported hardware

- g.USBamp (Windows only)
- g.Mobilab
- g.BSamp
- BrainProducts Brainamp series (Windows only)
- generic joysticks
- software sine generator
- LifeTool IntegraMouse
- Generic mice

## 1.2 Planned hardware support

- National Instruments DAQ cards
- Generic Keyboards
- NIRx NIRScout
- Neurosky MindSet
- Tunable EEG simulator

## 1.3 Planned extensions

- remote configuration
- client-based channel and sampling rate selection
- acquired data storage using .gdf files
- streaming of stored files

## 1.4 License

The TOBI signal server is licenced under the [GPLv3](#).

## 1.5 Contact

For further information please contact [SignalServer@tobi-project.org](mailto:SignalServer@tobi-project.org).

## 2 Hardware Requirements

- Run TiA with highest process priority possible PC requirements:
- CPU: at least 200 MHz (already tested on embedded systems)
- RAM: 32 MB

For networking usage:

- Ethernet min. 100 MBit (1 GBit recommended)  
(needed network connection varying by the sampling rate and the number of channels acquired)

## 3 Installation Instruction

### 3.1 Hardware Drivers

To use data acquisition devices as g.tec's g.USBamp, the respective drivers, provided by the manufacturer have to be installed first.

**Drivers provided by the manufacturers:**

- g.USBamp (Windows only)
- BrainProducts Brainamp series (Windows only)

### 3.2 Debian/Ubuntu

Download `libtia.deb` and `signalserver.deb` for your respective operating system and platform (32 or 64 bit). Install it with your preferred packet manager (e.g. synaptic). `libtia.deb` has to be installed first.

Done ;-)

*All required files are available at [www.tobi-project.org/download](http://www.tobi-project.org/download).*

### 3.3 Windows

Download the Signal Server setup file and install it.

Done ;-)

*All required files are available at [www.tobi-project.org/download](http://www.tobi-project.org/download).*

## 4 Using the Signal Server

### 4.1 Generic Usage

The Signal Server is shipped with two demo configuration files ("`server_config_comments.xml`" and "`server_config.xml`"). Those files are showing exemplary configurations of the Signal Server. The "`_comments.xml`" file is equipped with additional comments to facilitate understanding the configuration of the Signal Server.

The Signal Server can be started by a quickstart, where the file "`server_config.xml`" is automatically read (operating system dependent) or by passing an individual config file to the Signal Server.

**Start Commands for individual config files:**

```
signalserver your_config_file.xml
```

o r

```
signalserver -f your_config_file.xml
```

### **Commands while the Signal Server is running:**

q ... stop

r ... restart (problems with various DAQ devices – bugfix in progress)

## **4.2 Debian/Ubuntu**

With the first start of the Signal Server, a folder named “tobi\_sigserver\_cfg” is automatically created within the users home folder. Two files named “server\_config\_comments.xml” and “server\_config.xml” are located inside this folder.

Starting the Signal Server without any parameters automatically reads and starts the “server\_config.xml” inside the “tobi\_sigserver\_cfg” folder (quickstart).

## **4.3 Windows**

Starting the Signal Server without any parameters (quickstart) the program automatically reads and starts the “server\_config.xml” inside the folder where the Signal Server is installed (e.g. C:\Program Files\TOBI SignalServer\server\_config.xml).

To use a different configuration, an other .xml file can be handed to the Signal Server either using the commandline or just dragging the respective file on the Signal Server executable or a link to it.

# **5 The XML Configuration File – Basics**

## **5.1 Subject Information**

## **5.2 Server Settings**

## **5.3 The Hardware sections**

### **Mode**

- Master
- Slave
- Aperiodic

### **Device Settings**

- Sampling Rate
- Blocksize
- Measurement Channels

## **Channel Settings**

- Selection

## **6 Specific Hardware Configuration with the XML Config File**

### **6.1 Sine Generator**

### **6.2 g.USBamp**

### **6.3 g.Mobilab**

### **6.4 g.BSBamp**

### **6.5 Brainamp Series**

### **6.6 EEG Simulator**

### **6.7 Generic Joysticks**

### **6.8 IntegraMouse + Generic Mouses**

---